sifin diagnostics gmbh | Berliner Allee 317-321 | 13088 Berlin | Germany | T: +49 30 927 030-0 | F: +49 30 927 030-30 | info@sifin.de | www.sifin.de

Anti-Salmonella Gr. B, C, D, E
Anti-Salmonella O
Anti-Salmonella Vi

de

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Anti-Salmonella Gr. B, C, D, E
Anti-Salmonella O
Anti-Salmonella Vi

en

Test reagent for slide agglutination

INFORMATION FOR PROFESSIONAL USE

CE IVD

Intended use
The Anti-Salmonella test reagents are intended for use either in identifying group-specific O-antigens, and thus for the assignment of an isolate to the corresponding group of the Kauffmann-White scheme, or in the detection of additional O-antigens, insofar as these are required to determine or check the serovar. Additional serological differentiation must be performed using the Anti-Salmonella H test reagent in accordance with the White-Kauffmann-Le Minor scheme (Kauffmann-White scheme).

Principle of the test
If the *Salmonella* strain possesses an antigen that is covered by the test reagent, this antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

Composition
The test reagents are monoclonal antibodies, test sera or a mixture of monoclonal antibodies and test serum. Monoclonal antibodies are produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the respective Salmonella antigens.
The test reagents are immunized rabbit sera from which unspecific agglutinins have been removed by absorption.
Preservative: Sodium azide (NaN_3) 0.9 mg/ml

Available specificities Anti-Salmonella

Specificity	REF	Vol	Form
Gr. B	TR 1201	1 ml	liquid
	TR 5201	5 ml	liquid
Gr. C	TR 1202	1 ml	liquid
Gr. D	TR 1203	1 ml	liquid
	TR 5203	5 ml	liquid
Gr. E	TR 1204	1 ml	liquid
O:2	TR 1301	1 ml	liquid
O:4	TR 1302	1 ml	liquid
	TR 5302	5 ml	liquid
O:5	TR 1303	1 ml	liquid
	TR 5303	5 ml	liquid
O:6_1 [1]	TR 1304	1 ml	lyophilized
O:7	TR 1305	1 ml	liquid
O:8	TR 1306	1 ml	liquid
O:9	TR 1307	1 ml	liquid
	TR 5307	5 ml	liquid
O:10	TR 1308	1 ml	liquid
O:11	TR 1323	1 ml	lyophilized
O:13	TR 1325	1 ml	lyophilized
O:14 [2]	TR 1309	1 ml	lyophilized
O:15	TR 1310	1 ml	liquid
O:16	TR 1328	1 ml	lyophilized
O:17	TR 1329	1 ml	lyophilized
O:18	TS 1330	1 ml	lyophilized
O:19 [5]	TR 1311	1 ml	liquid
O:20	TR 1312	1 ml	liquid
O:21	TR 1331	1 ml	lyophilized
O:22	TS 1332	1 ml	lyophilized
O:25	TR 1335	1 ml	lyophilized
O:27	TR 1313	1 ml	liquid
O:28	TR 1336	1 ml	lyophilized
O:30	TR 1339	1 ml	lyophilized

Specificity	REF	Vol	Form
O:34 [3]	TR 1314	1 ml	liquid
O:35	TR 1341	1 ml	lyophilized
O:38	TR 1344	1 ml	lyophilized
O:39	TR 1345	1 ml	lyophilized
O:40	TR 1346	1 ml	lyophilized
O:41	TR 1347	1 ml	lyophilized
O:42	TR 1348	1 ml	lyophilized
O:43	TR 1349	1 ml	lyophilized
O:44	TR 1350	1 ml	lyophilized
O:45	TR 1351	1 ml	lyophilized
O:46	TR 1315	1 ml	liquid
O:47	TR 1353	1 ml	lyophilized
O:48	TR 1354	1 ml	lyophilized
O:50	TR 1355	1 ml	lyophilized
O:51	TR 1356	1 ml	lyophilized
O:52	TR 1357	1 ml	lyophilized
O:53	TR 1358	1 ml	lyophilized
O:54	TR 1359	1 ml	lyophilized
O:55 [4]	TR 1360	1 ml	lyophilized
O:56	TR 1361	1 ml	lyophilized
O:57	TR 1362	1 ml	lyophilized
O:58	TR 1363	1 ml	lyophilized
O:59	TR 1364	1 ml	lyophilized
O:60	TR 1365	1 ml	lyophilized
O:61	TR 1366	1 ml	lyophilized
O:62	TR 1367	1 ml	lyophilized
O:63	TR 1368	1 ml	lyophilized
O:65	TR 1369	1 ml	lyophilized
O:66	TR 1370	1 ml	lyophilized
O:67	TR 1371	1 ml	lyophilized
Vi	TR 1316	1 ml	liquid

[1] To be used for testing serovars of the group C_2-C_3 for the presence or absence of the factor O:6.
[2] Covers O:14 of the group H but not of the group C_2.
[3] May only be used for differentiation within group E, since it crossreacts with O:12_2, O:28 and O:55.
[4] Reacts with some strains of the group O:59.
[5] Reacts with some strains of the group B with O:27


Form in which product is supplied, shelf life and conditions for storage
Test reagents that are frequently needed in the microbiological laboratory are available in liquid form ready to use. Test reagents for specificities that are needed less often are lyophilized (freeze-dried). If stored at 2...8 °C both before and after opening, **liquid** test reagents may be used up to the date given on the label. The reagents are ready to use. After use, the bottle must be properly closed.
Lyophilized test reagents must be dissolved in 1 ml distilled water before use. They must be properly closed with the enclosed pipette. If stored at 2...8 °C, they remain usable for at least 18 months. However, they must not be used after the date given on the label. The test reagent may sometimes show turbidity not caused by microbes. Such turbidity does not impair effectivity and the test reagents can be clarified by centrifugation or filtration.
The test reagents must have their temperatures adjusted to room temperature (18...26 °C) before use.

Warnings and precautions
The biotechnological manufacture of the monoclonal antibodies means that the risk of contamination by infectious agents can be virtually excluded. Because they contain animal materials (fetal calf serum, stabilizer), they should be treated as potentially infectious and handled accordingly. Test reagents containing biological material in the form of rabbit serum should therefore be treated as potentially infectious and handled accordingly. As these materials contain sodium azide, contact with the skin and mucous membranes must be avoided. In case of contact, rinse with plenty of water.
Since the performance of the slide agglutination test involves working with native pathogenic materials, all necessary work protection procedures must be adhered to (risk of infection)!


Materials and equipment not supplied
Glass slides, stirring rods, physiological saline (NaCl 9 g/l), disposal containers for infectious material, White-Kauffmann-Le-Minor scheme.

Test material and methodology
Transfer a small amount of bacterial mass from a 16-20 hour old subculture (e.g. nutrient agar, blood agar or Kligler culture medium) onto a slide and mix well with a drop of test reagent (ca. 25 µl), so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface. The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth). To control for spontaneous agglutination, a negative control (NC) must be performed at the same time using physiological saline in place of the test reagent.
Note: If the Vi-antigen is present, it can mask the identity of O-antigens. To detect the O-antigens, it may be necessary to heat an antigen suspension of the isolate for 60 minutes at 100 °C or for 15 minutes at 120 °C.

Evaluation
The test can only be evaluated if the negative control (NC) remains milky-opaque.
Positive: visible agglutination after the sample has been tilted back and forth less than 20 times. In a strongly positive reaction (PR), agglutination (coarsely or finely flocculent) appears as soon as the bacterial mass is mixed in. In a weakly positive result, agglutination only appears after the slide has been tilted back and forth 10-20 times.

Negative: if the suspension remains milky-opaque or the reaction begins to occur only after the slide has been tilted back and forth more than 20 times, the result is negative (NR).

Quality assurance during the testing procedure
For quality control of serotyping using slide agglutination, it is important that the strains used for this purpose express their antigens well on the cell surface. It is therefore recommended to use strains from interlaboratory tests, field strains of defined origin that have been fully characterised in an external laboratory or suitable commercial *Salmonella* O- and H-antigens.

Limitations of the method
The test reagents react solely with *Salmonella* strains which contain antigens of the specificity named in the declaration. In exceptional cases, there is the possibility of crossreactions with other genera of the Enterobacteriaceae family (e. g. with strains of *Citrobacter* spp., *Hafnia alvei, Proteus* spp. or *E. coli*) due to antigen identities or related antigens. For this reason, biochemical examination is required to establish membership of the genus Salmonella. For further serological differentiation, the Anti-Salmonella H test reagent is used according to the White-Kauffmann-Le Minor scheme.

Explanation of the symbols used

LOT Batch code (Lot)
REF Catalogue number
IVD In Vitro Diagnostic Medical Device
TR Test reagent
mTR Monoclonal test reagent
Use by YYYY-MM (MM = end of month)
Temperature limitation
Consult instructions for use
Slide agglutination
LYO lyophilized

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