sifin diagnostics gmbh | Berliner Allee 317-321 | 13088 Berlin | Germany | T: +49 30 700 144-0 | F: +49 30 700 144-30 | info@sifin.de | www.sifin.de

Oxidase-Reagenz

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Oxidase-Reagenz	mg
N,N,N',N'-Tetramethyl-1,4-phenylendiamindihydrochlorid	20,0

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en

REF TN1517 Vol. 12 x 1 bottle

Intended use
Used to test for the oxidase reaction.

Composition
1 bottle of Oxidase Reagent contains:

Oxidase Reagent	mg
Potassium N,N,N′,N′-Tetramethyl-p-phenylenediamine dihydrochloride	20.0

Shelf life and storage
If stored at 15...25 °C, Oxidase Reagent may be used up to the date given on the label. It must be stored protected from light.

Warnings and precautions
Classification according to Regulation (EC) 1272/2008 (CLP):

Signal word: Warning

Hazard elements:
Contains Potassium N,N,N′,N′-Tetramethyl-p-phenylenediamine dihydrochloride, CAS 637-01-4

H-statements: H315: Causes skin irritation. H319: Causes serious eye irritation. H335: May cause respiratory irritation.
P-statements: P261 - Avoid breathing dust, spray. P302+P352 - IF ON SKIN: Wash with plenty of water. P304+P340 - IF INHALED: Remove person to fresh air and keep comfortable for breathing. P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P332+P313 - If skin irritation occurs: Get medical advice/attention. P337+P313 - If eye irritation persists: Get medical advice/attention.

Directions
Dissolve the contents of one bottle Oxidase Reagent in 2 ml sterile distilled water.

Instrucions for use
Put 2-3 drops onto a piece of filter paper in a Petri dish. Triturate some colony material on the prepared filter paper using a platinum, plastic (not nickel-chrome) wooden or glass rod. The appearance of a deep blue colour within 30 seconds is counted as a positive reaction.

Explanation of the symbols used

Batch code (Lot)

Catalogue number

Consult instructions for use

Use by YYYY-MM (MM = end of month)

Temperature limitation

Date of revision: 29/04/2025