sifin diagnostics gmbh | Berliner Allee 317-321 | 13088 Berlin | Germany | T: +49 30 927 030-0 | F: +49 30 927 030-30 | info@sifin.de | www.sifin.de

Anti-Salmonella Poly-H Phase 1 & 2

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Anti-Salmonella Poly-H Phase 1 & 2

Monoclonal test reagent for slide agglutination

INFORMATION FOR PROFESSIONAL USE

CE IVD

en

	REF	Vol.
Anti-Salmonella Poly-H Phase 1 & 2	TR 1141	1 ml
	TR 1145	5 ml

Intended use
The test reagent is intended for use in the serological identification of the H-antigens of *Salmonella* strains isolated from test material of varied origin, using slide agglutination.

Principle of the test
If the isolated strain contains a *Salmonella*-H-antigen that is covered by the test reagent, this antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

Composition
The Anti-Salmonella Poly-H Phase 1 & 2 test reagent is a mixture of monoclonal antibodies. It is produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the respective group specific H-antigens.
The Anti-Salmonella Poly-H Phase 1 & 2 test reagent contains antibodies for the detection of the following H-antigens or H-antigen complexes: a, b, c, d, E, G, i, k, L, r, y, z, Z_4, Z_6, Z_{10}, Z_{29}, Z_{35}, Z_{38}, Z_{41} and H:1. Antibodies which detect complex antigens, e.g. Anti-H:E, H:G, H:L or H:1, detect all combinations of the antigens H:e, H:g (including H:m,t), H:ℓ, or H:1 respectively.
Preservative: Sodium azide (NaN_3) 0.9 mg/ml

Form in which product is supplied, shelf life and conditions for storage
The Anti-Salmonella Poly-H Phase 1 & 2 test reagent is available in liquid form ready to use.
If stored at 2...8 °C before and after opening, it may be used up to the date given on the label. The reagent is ready to use. After use, the bottle must be properly closed.
In a few rare cases, slight turbidity not caused by microbes may occur during storage. Such turbidity does not impair effectivity and the test reagents can be clarified by centrifugation or filtration.
The test reagents must have their temperatures adjusted to room temperature (18...26 °C) before use.

Warnings and precautions
The biotechnological manufacture of the test reagent means that the risk of contamination by infectious agents can be virtually excluded. Because it contains animal materials (fetal calf serum, stabilizer), it should be treated as potentially infectious and handled accordingly.
The test serum is a biological material which should be treated as potentially infectious and handled accordingly.
As this material contains sodium azide, contact with the skin and mucous membranes must be avoided. In case of contact, rinse with plenty of water. Since the performance of the slide agglutination test involves working with native pathogenic materials, all necessary work protection procedures must be adhered to (risk of infection)!


Materials and equipment not supplied
Black worktop, suitable light source, glass slides, stirring rods, physiological saline (NaCl 9 g/l), sterile distilled water, swarm agar, sterile petri dishes (6 cm diameter), disposal containers for infectious material, White-Kauffmann-Le-Minor scheme (Kauffmann-White scheme).

Test material and methodology
1. Take a small amount of bacterial mass from a 16-20 hour old sub-culture (e.g. Kligler culture medium), picked from the lower, moister layer of the agar slope. Transfer onto a slide and mix well with one drop (ca. 25 µl) of Anti-Salmonella Poly-H Phase 1 & 2 test reagent, so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface.
The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth). To control for spontaneous agglutination, a negative control must be performed at the same time using physiological saline in place of the test reagent.
2. In case the H-antigen test performed according to the method described above is not successful, the strain must be inoculated onto swarm agar, to enhance the development of the flagella antigens. For the swarm agar we recommend our ready-to-use culture medium sifin REF TN 1702.
After liquefaction, place 10 ml of this swarm agar into a petri dish (6 cm diameter). Once it has solidified, sprinkle the surface of the agar with ca. 100 µl sterile distilled water and apply the strain in the form of a point in the centre of the plate. Incubate the dish overnight with the agar layer facing down at 35...37 °C. For the slide agglutination test, collect material from the edge of the dish.

Evaluation
The test can only be evaluated if the negative control (NC) remains milky-opaque.
Positive: visible agglutination after the sample has been tilted back and forth less than 20 times. In a strongly positive reaction (PR), agglutination (coarsely or finely flocculent) appears as soon as the bacterial mass is mixed in. In a weakly positive result, agglutination only appears after the slide has been tilted back and forth 10-20 times.

Negative: if the suspension remains milky-opaque or the reaction begins to occur only after the slide has been tilted back and forth more than 20 times, the result is negative (NR).

In case of a positive reaction, biochemical confirmation of the genus *Salmonella* must be obtained before further serological differentiation of the strain can be carried out using poly-specific, group-specific and monospecific Anti-Salmonella test reagents in accordance with the White-Kauffmann-Le-Minor scheme. The only non-flagellated *Salmonella* serovar in the White-Kauffmann-Le-Minor scheme (*S*. Gallinarum) cannot be detected using the Anti-Salmonella H test reagents.

Quality assurance during the testing procedure
For the quality control of serological determination of *Salmonella* H antigens by slide agglutination test, the good expression of the strain's cell surface antigen is important. The use of strains from interlaboratory tests, or field strains of defined origin that have been fully characterised by an external laboratory, or the use of suitable commercially available *Salmonella* test antigens is therefore recommended for quality control.

Limitations of the method
The test reagent reacts solely with *Salmonella* strains which contain antigens of the specificity named in the declaration. Cross-reactions with strains of other genera of *Enterobacteriaceae* have not been found.
On its own, positive serological results obtained using Anti-Salmonella Poly-H test reagent do not constitute proof of the presence of Salmonella. For this reason, biochemical means must be used to confirm that the strain belongs to the *Salmonella* genus. A negative result does not exclude the presence of Salmonella with absolute certainty.

Explanation of the symbols used

LOT	Batch code (Lot)
REF	Catalogue number
IVD	In Vitro Diagnostic Medical Device
mTR	Monoclonal test reagent
	Use by YYYY-MM (MM = end of month)
	Temperature limitation
	Consult instructions for use
	Slide agglutination

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