sifin diagnostics gmbh | Berliner Allee 317-321 | 13088 Berlin | Germany | T: +49 30 700 144-0 | F: +49 30 700 144-30 | info@sifin.de | www.sifin.de

Anti-D

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Anti-D

Monoclonal test reagent for blood-group determination

INFORMATION FOR PROFESSIONAL USE

CE 0483 IVD

en

Intended use

Anti-D is an in vitro diagnostic and is used for blood grouping. It is employed by specialized personnel to qualitatively detect RhD antigens of the rhesus blood type system by erythrocyte agglutination in human blood samples [with or without anticoagulants (EDTA or PAGGS-M)] in manual plate and tube tests. The category VI partial RhD antigen is not detected with this clone.

Principle of the test

Along with the antigens of the ABO blood group system, antigen D is the most important blood group antigen in blood group serology. Unlike those of the ABO system, antibodies of the Rh system do not occur generally in the serum. They are always the result of exposure to the antigen during pregnancy or through transfusion.
The Rh characteristic D occurs in different forms of expression:

- as a completely expressed characteristic, typical for Rh positive (D positive)
- as a characteristic that is weakened and/or changed in expression: either as weakly expressed antigen D (weak D) or as qualitatively changed antigen D (partial D) which can also be weakly expressed (e. g. D category VI = D^{VI}). Persons with a significantly qualitatively changed D antigen, especially persons with the category D^{VI}, can be immunized by a completely expressed Rh characteristic D.

The guidelines require that the testing of patients, pregnant women and newborns be performed using at least two Anti-D test reagents with different IgM antibodies which do not agglutinate the category D^{VI}.
The testing of blood donors must also be performed using two Anti-D test reagents. It is necessary that every weakly or only partially expressed Rh characteristic antigen D is reliably detected (see „Interpretation of test results“).
The procedures used with the test reagent Anti-D are based on the principle of agglutination. Agglutination of the red blood cells is a positive test result and shows the presence of the Rh antigen D. The presence of no agglutination constitutes a negative result and indicates that Rh antigen D is not detectable.

Composition

The test reagent is manufactured from the cell culture supernatant of the heterohybridoma cell line BS225. The antibodies belong to the human IgM isotype. The addition of stabilizing substances ensures that the test reagent remains active throughout its shelf life. Preservative: sodium azide (≤0.99 mg/ml).

Shelf life and storage

If stored at 2...8 °C, the test reagent can be used until the date indicated on the label. After the first opening, it must be stored well closed at 2...8 °C. Avoid bacterial contamination. Do not use if markedly turbid.

Test material

Blood should be drawn by an aseptically, with or without an anticoagulant (EDTA or PAGGS-M). The specimen should be tested as soon as possible after collection. If testing is delayed, the specimen must be stored at 2...8 °C.
Bacterial contamination of the specimen may cause false test results. Blood drawn into EDTA as well as clotted blood may be tested up to 14 days after collection, and donor blood may be tested up to the expiration date. Note, however, that storage may result in weaker than normal reactions (see „Error sources and limitations“ below).

Materials and equipment not supplied

Test tubes (75 × 12 mm or 70 × 8 mm), plates, pipettes (drop volume: 40-50 µl), applicator sticks, physiological sodium chloride solution (154 mmol/l), laboratory centrifuge, timer.

Methods

1. Plate test

1.1. Place 1 drop (40-50 µl) of test reagent Anti-D on a properly labelled plate.
1.2. To each drop of test reagent, add 1 drop (40-50 µl) of a 5 % suspension of the red blood cells from the blood to be tested (in physiological sodium chloride solution). Mix carefully with a clean glass or plastic rod.
1.3. Read for agglutination after 5, at the most 20 min incubation at 18...26 °C by rotating or rocking the plate. Alternatively, incubate for 20 min at 37 °C.
1.4. Record test results. Test results must be interpreted immediately upon completion of the test as the reaction product is unstable.

2. Tube test

2.1. Prepare a 3-5 % suspension of the red blood cells to be tested (washed or unwashed) in physiological sodium chloride solution.
2.2. Place 1 drop (40-50 µl) of the test reagent Anti-D in a small, properly labelled test tube (e. g., 70 × 8 mm).
2.3. Add 1 drop (40-50 µl) of the suspension to each tube and mix thoroughly by shaking.
2.4. Centrifuge immediately for 1 to 2 min at rcf 180 to 270 × *g*.
2.5. Gently shake the tube to loosen the sediment from the bottom of the tube and read for agglutination.
2.6. Record test results. Test results must be interpreted immediately upon completion of the test as the reaction product is unstable.
2.7. In the case of a weak positive, doubtful or negative result, homogenize the tube contents by shaking; incubate the tube for 15 to 20 min at 18...26 °C and centrifuge and read again (see 2.4. to 2.6.)

Interpretation of the test results

Interpretation of the test results of patients, pregnant women and newborns: In the case of concurring positive reactions of the two anti-D test reagents, including unambiguously weak positive reactions, the result is “Rh positive (D positive)”. If the two test reagents yield different or questionable positive results, the patient is to be declared as „recipient Rh negative (D negative)“.
Interpretation of the results of blood donors: In the case of concurring positive reactions of the two anti-D test reagents, including unambiguously weak positive reactions, the result is „Rh positive (D positive)“
In all other cases the sample must be examined using suitable polyclonal or oligoclonal reagents in the indirect anti-human globulin test or a suitable RhD PCR must be carried out in order to also be able to detect any weak or only partial D antigen (see guidelines listed above or country-specific regulations).

Warnings and precautions

- The use of biotechnological production means that the risk of contamination of the test reagent with infectious agents is virtually excluded. Since the test reagent contains materials of animal origin (fetal calf serum, stabilizer), it should be treated as a potential health hazard and be handled accordingly.
- The test reagent contains sodium azide; avoid any contact with the skin or mucous membranes.
- Materials from processed tests must be disposed of in accordance with the appropriate national guidelines, e.g. in Germany as waste of Waste Catalog Ordinance (AVV) group 18 01 (AS 18 01 03) as set out in the “Guideline on the proper disposal of waste from health-care establishments“ (LAGA guideline, current version).

Notes and limits of the procedure

The performance evaluation according to the Commission Implementing Regulation (EU) 2022/1107 laying down common specifications (CS) yielded the following calculated diagnostic values:

	Σ	positive	false negative	sensitivity	negative	false negative	specificity
Tube test	1430 [1]	1142	0	100 % [1]	288	0	100 %
	1516 [2]	1221 [2]	7 [2]	99.43 % [2]	292 [2]	0	100 % [2]
Plate test	2168 [1]	1686	1	99.94 % [1]	524	0	100 %
	2250 [2]	1726 [2]	39 [2]	97.79 % [2]	527 [2]	0	100 % [2]

1) Samples with normal antigen expression 2) including weak D samples

- Weakly expressed D antigen (weak D) is recognized by the plate test with only insufficient sensitivity. Therefore, the tube test with post-incubation in the case of negative, doubtful or weak positive results is recommended.
- A weakly positive result may indicate the presence of a weak D or a partial D (of one category).
- Red blood cells which do not possess the D antigen are not agglutinated.
- Of the more frequent categories D^{II} to D^{VII}, the only one with which the monoclonal IgM antibody BS225 does not react is the category D^{VI}.

Error sources and limitations

- **Controls:** For each test, positive and negative control red blood cells and an auto control (test for autoagglutination) must be tested in parallel.
- Older blood specimens may show weaker reactions than fresh specimens.
- The concentrations of red blood cell suspensions and the reaction times and conditions specified above must be observed in order to obtain correct results.
- Inappropriate storage decreases the activity of the reagent.
- Titre determination: To prepare dilutions of the monoclonal test reagents the use of physiological saline containing 2 % protein (bovine serum albumin) is recommended to facilitate the resuspension of the reaction products from the plate or the bottom of the tube.
- The application procedures are for the described manual methods only. When using automated or semi-automated instruments, laboratories must follow the instructions contained in the operator‘s manual provided by the manufacturer of the instrument and perform validations according to approved procedures.

Vigilance

All serious incidents occurring in connection with the product must be reported to sifin diagnostics gmbh as manufacturer and the responsible authority of the member state in which the user and/or the patient is based. The current summary report on safety and performance is available on request (info@sifin.de) or via EUDAMED.

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Explanation of the symbols used

LOT Batch code (Lot)
REF Catalogue number
IVD In Vitro Diagnostic Medical Device
Rh Blood grouping test reagent
ACT Activity/ Minimum titre
CLONE Clone designation

Use by YYYY-MM (MM = end of month)
Temperature limitation
www.sifin.de Consult instructions for use

Additional languages of the instruction for use of your product you will find here: www.sifin.de

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