sifin diagnostics gmbh | Berliner Allee 317-321 | 13088 Berlin | Germany | T: +49 30 927 030-0 | F: +49 30 927 030-30 | info@sifin.de | www.sifin.de

# Anti-Salmonella H

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# Anti-Salmonella H

Test reagent for slide agglutination and phase reversal (PR)

INFORMATION FOR PROFESSIONAL USE

CE IVD

en

**Intended use**
The Anti-Salmonella H test reagents are intended for use in identifying or checking the H-antigens of *Salmonella* strains in accordance with the White-Kauffmann-Le-Minor scheme (Kauffmann-White scheme) using slide agglutination. They enable the serovar to be determined.
The test reagents indicated in the table below may also be used in the case of diphasic strains for induction of the unidentified H-phase (phase inversion by culture on swarm agar/ Sven Gard method).

**Principle of the test**
If the *Salmonella* strain contains an H-antigen that is covered by the test reagent, this antigen becomes bound when mixed with the specific antibody. The antigen-antibody reaction results in clearly visible agglutination of the strain.

**Composition**
The test reagents are monoclonal antibodies, test sera or a mixture of monoclonal antibodies and test serum. Monoclonal antibodies are produced from cell culture supernatants of hybridoma cell lines which secrete antibodies against the respective Salmonella H antigens. The test sera are immunized rabbit sera from which unspecific agglutinins have been removed by absorption.
Preservative: Sodium azide ($NaN_3$) 0.9 mg/ml

**Available specificities Anti-Salmonella**

| Specif. | Antigens identified | PR |
|---|---|---|
| H:a | H:a | x |
| H:b | H:b | x |
| H:c | H:c | x |
| H:d | H:d | x |
| H:E | H:e,h; e,n,x; e,n,z15; e,n,x,z15 | x |
| H:f | H:f,g; f,g,s; f,g,t | |
| H:g | All antigen combinations of H:G complexes, except H:m,t | x |
| H:g,m | All antigen combinations of H:G complexes, incl. H:m,t | |
| H:h | H:e,h | |
| H:i | H:i | x |
| H:k | H:k | x |
| H:L | H:ℓ,v; ℓ,w; ℓ,z13; ℓ,z28; ℓ,z13,z28 | x |
| H:m | H:g,m; g,m,s; g,m,s,t; g,m,q; g,m,p,s; g,m,t; m,p,t,u; m,t | |
| H:n | H:e,n,x; e,n,z15; e,n,x,z15 | x |
| H:p | H:g,m,p,s; g,p; g,p,s; g,p,u; m,p,t,u | |
| H:q | H:g,q; g,m,q | |
| H:r | H:r | x |
| H:s | H:f,g,s; g,m,s; g,m,s,t; g,p,s; g,s,t | |
| H:t | H:m,t; f,g,t; g,m,t; g,m,s,t; g,t; g,s,t; m,p,t,u | |
| H:u | H:g,p,u; m,p,t,u | |
| H:v | H:l,v | |
| H:w | H:l,w | |
| H:x | H:e,n,x; e,n,x,z15 | |
| H:y | H:y | x |
| H:z [1] | H:z (I, II, III) | x |
| H:z4,z23 | H:z4,z23; z4,z24; z4,z32; z4,z23,z32 | |
| H:z6 | H:z6 | |
| H:z10 | H:z10 | x |
| H:z15 | H:e,n,z15; e,n,x,z15 | |
| H:z23 | H:z4,z23; z4,z23,z32 | |
| H:z24 | H:z4,z24 | |
| H:z28 | H:ℓ,z28; ℓ,z13,z28 | |
| H:z29 | H:z29 | |
| H:z32 | H:z4,z32 | |
| H:z35 | H:z35 | x |
| H:z38 | H:z38 | x |
| H:z41 | H:z41 | x |
| H:1 | H:1,2; 1,5; 1,6; 1,7; 1,2,7; 1,5,7 | x |
| H:2 | H:1,2 | |
| H:5 | H:1,5 | |
| H:6 | H:1,6 | |
| H:7 | H:1,7 | |

[1] identifies Hz in Subspecies I, II und III

**Form in which product is supplied, shelf life and conditions for storage**
Test reagents that are frequently needed in the microbiological laboratory are available in liquid form ready to use. Test reagents for specificities that are needed less often are lyophilized (freeze-dried).
If stored at 2...8 °C both before and after opening, **liquid** test reagents may be used up to the date given on the label and must be properly resealed after use. The lyophilized test reagents must be dissolved in 1 ml or 5ml distilled water in accordance with the declaration before use. They can be properly closed with the enclosed pipette. If stored at 2...8 °C, they remain usable for at least 18 months. However, they must not be used after the date given on the label.
The test reagent may sometimes show turbidity not caused by microbes. Such turbidity does not impair effectivity and the test sera/ test reagents can be clarified by centrifugation or filtration.
Prior to use, the test reagents must be warmed up to room temperature (18...26 °C).

**Warnings and precautions**
The biotechnological manufacture of the test reagents means that the risk of contamination by infectious agents can be virtually excluded. Because they contain animal materials (fetal calf serum, stabilizer), they should be treated as potentially infectious and handled accordingly. Test reagents containing biological material in the form of rabbit serum should therefore be treated as potentially infectious and handled accordingly.
As these materials contain sodium azide, contact with the skin and mucous membranes must be avoided. In case of contact, rinse with plenty of water. Since the performance of the slide agglutination test involves working with native pathogenic materials, all necessary work protection procedures must be adhered to (risk of infection)!

**Materials and equipment not supplied**
Glass slides, stirring rods, physiological saline (NaCl 9 g/l), disposal containers for infectious material, Kauffmann-White-Le-Minor scheme, sterile petri dishes (6 cm diameter), swarm agar, microwave oven or water-bath, thermometer.

**Test material and methodology**
**1. Determination of the antigen**
**1.1** Take a small amount of bacterial mass from a 16-20 hour old subculture (e.g. nutrient agar or Kligler culture medium), picked from the lower, moister layer of the agar slope. Transfer onto a slide and mix well with a drop (ca. 25 µl) of Anti-Salmonella H test reagent, so that a homogenous, slightly milky suspension results. The slide should be placed on a dark surface.
The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth). To control for spontaneous agglutination, a negative control must be performed at the same time using physiological saline in place of the test reagent.
**1.2** In case the H-antigen test performed according to the method described above is not successful, the strain must be inoculated onto swarm agar, to enhance the development of the flagella antigens. For the swarm agar we recommend our ready-to-use culture medium sifin REF TN 1702. After liquefaction, place 10 ml of this swarm agar into a petri dish (6 cm diameter). Once it has solidified, sprinkle the surface of the agar with ca. 100 µl sterile distilled water and apply the strain in the form of a point in the centre of the plate. Incubate the dish overnight with the agar layer facing down at 35...37 °C. For the slide agglutination test, collect material from the edge of the dish.
**2. Phase reversal**
If the strain is diphasic but only one H phase can be identified, the phase inversion procedure (Sven Gard method) must be performed to induce the second phase.
**2.1.** Mix 0.1-0.2 ml Anti-Salmonella H test reagent with 10 ml liquefied swarm agar (cooled to 40...45 °C) in a 6 cm diameter petri dish. After the agar has solidified, sprinkle the plate with distilled water (as described under 1.2.), inoculate it with the Salmonella strain, incubate it and perform the slide agglutination test. If the phase inversion is unsuccessful, the procedure should be repeated.
**2.2.** As an alternative to this classical procedure, we recommend the following simplified procedure. Add 0.1-0.2 ml Anti-Salmonella H test reagent drop by drop to a plate (6 cm diameter) containing 10 ml swarm agar, and spread over the surface with a sterile glass spatula. Then inoculate the middle of the plate and incubate as indicated above. After 16-20 hours, the strain which swarms can be used to identify the second H phase. This simplified procedure usually enables the second phase to be identified at the first attempt.

**Evaluation**
The test can only be evaluated if the negative control remains milky-opaque.
**Positive:** visible agglutination after the sample has been tilted back and forth less than 20 times. In a strongly positive reaction, agglutination (coarsely or finely flocculent) appears as soon as the bacterial mass is mixed in. In a weakly positive result, agglutination only appears after the slide has been tilted back and forth 10-20 times.
**Negative:** if the suspension remains milky-opaque or the reaction begins to occur only after the slide has been tilted back and forth more than 20 times, the result is negative.

**Quality assurance during the testing procedure**
For quality control of the serotyping using slide agglutination, it is important that the strains used for this purpose express their antigens well on the cell surface. It is therefore recommended to use strains from interlaboratory tests, field strains of defined origin that have been fully characterised in an external laboratory or suitable commercial *Salmonella* O- and H-antigens.

**Limitations of the method**
The test reagents react solely with *Salmonella* strains which contain antigens of the specificity named in the declaration. Cross-reactions with strains of other genera of *Enterobacteriaceae* have not been found. The identification of the H-antigens is achieved by the serological identification of the O-antigen(s) and the biochemical allocation to the *Salmonella* genus.

**Explanation of the symbols used**

LOT Batch code (Lot)
REF Catalogue number
IVD In Vitro Diagnostic Medical Device
TR Test reagent
mTR Monoclonal test reagent
PI Phase reversal
Use by YYYY-MM (MM = end of month)
Temperature limitation
Consult instructions for use
Slide agglutination
LYO lyophilized

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