sifin diagnostics gmbh | Berliner Allee 317-321 | 13088 Berlin | Germany | T: +49 30 700 144-0 | F: +49 30 700 144-30 | info@sifin.de | www.sifin.de

Kontrollantigene Shigella

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CE IVD

en

Intended use
The control antigens are used to check the agglutinability of the sifin Anti-Shigella test reagents and for quality control when carrying out the slide agglutination test.

Principle of the test
Agglutination using an antigen-antigen-antibody reaction between specific *Shigella* antibodies and control antigens of characterised *Shigella* serovars.

Composition
The Shigella control antigens are killed bacterial suspensions of characterised *Shigella* serovars in phosphate-buffered saline (PBS) with 0.5 % formalin.

REF	Available specificities	Antigen mixture
TS 1510	Control antigen Shigella flexneri	*S. flexneri* Type 1b, 2a, 3a, 4a
TS 1511	Control antigen Shigella dysenteriae	*S. dysenteriae* Type 1, 2, 3, 7
TS 1512	Control antigen Shigella boydii	*S. boydii* Type 1, 2, 5, 8
TS 1513	Control antigen Shigella sonnei	*S. sonnei* S-form (phase I), *S. sonnei* F-form (phase II)

The package size is 2 ml.

Form in which product is supplied, shelf-life and storage
The control antigens are provided in liquid, ready-to-use form. If stored at 2...8 °C both before and after opening, they may be used up to the date given on the label. After use, the bottles must be properly closed.

Warnings and precautions
As these materials contain formalin, contact with the skin and mucous membranes must be avoided! In case of contact, rinse with plenty of water.

Materials and equipment not supplied
Glass slides, stirring rods, physiological saline (NaCl 9,0 g/l), disposal containers for infectious material, sifin test reagents

Test material and methodology
The control antigens are used in quality control of polyspecific sifin test Anti-Shigella reagents using the slide agglutination test.
The control antigen Shigella dysenteriae can also be used to test the monospecific test reagents Anti-Shigella dysenteriae type 1 and Anti-Shigella dysenteriae type 2.
The control antigen Shigella sonnei is also used to test the test reagent Anti-Shigella sonnei S-form (phase I), Anti-Shigella sonnei F-form (Phase II) and Anti-Shigella sonnei S-form & F-form (phases I & II).

Fields of application for quality control

Control antigen Shigella flexneri	TS1510	Control antigen Shigella dysenteriae	TS1511	Control antigen Shigella boydii	TS1512	Control antigen Shigella sonnei	TS1513
Anti-Shigella I	TR1811 TR1811-01	Anti-Shigella II	TS1821 TS1821-01	Anti-Shigella III	TS1831	Anti-Shigella I	TR1811 TR1811-01
Anti-Shigella flexneri	TS1901	Anti-Shigella dysenteriae Type 1	TS2001			Anti-Shigella sonnei S-form (phase I)	TR2012
		Anti-Shigella dysenteriae Type 2	TS2002			Anti-Shigella sonnei F-form (phase II)	TR2013
						Anti-Shigella sonnei S-& F-form (phases I&II)	TR2014

The control antigens must be carefully shaken prior to use until they are fully mixed.
One drop of each control antigen (about 25 µl) is mixed on a slide with one drop of the sifin test reagent (about 25 µl) in equal measure, as is required to conduct the serological test on the control antigen.
The slide should be placed on a dark surface. The reaction is read with the naked eye by holding the slide in front of a light source against a black background and rocking it (tilting it back and forth). To control for spontaneous agglutination, a negative control must be performed at the same time using physiological saline in place of the test reagent.

Evaluation of the agglutination result
Visible agglutination after the sample has been tilted back and forth less than 20 times. The negative control must still show a milky turbidity.
Note for the user: The control antigens are used in accordance with the above intended use exclusively for quality control of the serological typing with the slide agglutination test when using sifin test reagents.

Explanation of the symbols used

LOT Batch code (Lot)
REF Catalogue number
IVD In Vitro Diagnostic Medical Device
Slide agglutination
Use by YYYY-MM (MM = end of month)
Temperature limitation
Consult instructions for use

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